

# ddRAura™ S

---



**Value and performance**  
*without compromise*

**Swissray** 

# A multifunctional system *for more versatility*

The ddRAura™ S performs direct digital radiographic examinations on upright, seated and supine patients with unprecedented efficiency and image quality. The system includes motorized vertical arm and SID movement for fast and precise system set-up.




- ddRAura™ S features**
- Multifunctional system for all radiography examinations
  - Motorized height and SID adjustment
  - Effortless manual tube-arm, detector and tube rotation
  - Exceptional low X-ray beam for standing exams
  - Various generator power incl. energy assist
  - Robust design, maintenance friendly
  - Fits into very small examination rooms

ddRAura™ S includes a full field 43x43 cm flat panel detector and can be configured with different generator power incl. energy assist to be operated from a standard wall socket.

# A user interface *with a difference*

The unique ddRAura™ user interface simplifies radiographic procedures by automatically setting the required technical factors and parameters for the specific body part, view position and age group. The interface is designed to complete the examination setup with the fewest steps possible.



- 1 Patient Registration
- 2 Organ Selection
- 3 Procedure Mapping
- 4 Image Acquisition
- 5 Image Review and Transfer

The ddRAura™ S is the diagnostic tool of choice for any pulmonary disease such as tuberculosis and performs all kind of regular x-ray images. It includes the proven Swisray workstation, which provides fast and logical examination programming and system set-up for utmost positioning convenience and patient throughput.

Distributed by

Swissray Medical AG reserves the right to make changes in specifications and/or to discontinue any product at any time without prior notice or obligation and will not be liable for any consequences resulting from the use of this publication. Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

© 2017 Swissray Medical AG

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

SHM\_49\_700\_007\_01\_EN

**Swissray Medical AG**

Turbistrasse 25–27  
6280 Hochdorf  
Switzerland

Phone + 41 41 914 12 12  
Fax + 41 41 914 12 13  
sales@swissray.com  
www.swissray.com

**Swissray International, Inc.**

31 Gordon Road  
Piscataway, NJ08854  
USA

Phone + 1 800 903 5543  
Fax + 1 908 353 1237  
info@swissray.com  
www.swissray.com

**Swissray Asia Healthcare Co., Ltd.**

7F, N° 87 Xinhua 1st Rd.  
Taipei City 114  
Taiwan (R.O.C.)

Phone + 886-2-8792-2699  
Fax + 886-2-8792-1898  
info.asia@swissray.com  
www.swissray.com

**Swissray** 